

UNIVERSAL,
SIMPLE,
VERSATILE

ZD ZEST DENTAL SOLUTIONS
ZEST | DANVILLE MATERIALS | PERIOSCOPY

LOCATOR®

 **PRO-SCAN**

Belgium
Eikenenweg 71
3520 Zonhoven
Tel. + 32(0)11 822 650

The Netherlands
Postbus 90
6170 AB STEIN
Tel. +31(0)30 609 0095

Germany
Tel. +32(0)11 822 650

info@proscan.be
www.proscan.be

webstore.proscan.be

LOCATOR,
UNIVERSAL &
SIMPLE



AVAILABLE IN THE WEBSHOP:
webstore.proscan.be



- 1 - Locator in post & core structure
- 2 - Locator with frame
- 3 - Locator on bar constructions
- 4 - Customized, individual Locator via CAD/CAM



distributed by
 **PRO-SCAN**



Self-aligning design


LOCATOR IMPLANT SYSTEM

Available for multiple implants systems, in multiple gingival heights



LOCATOR ROOT SYSTEM

For placement in a root canal



LOCATOR BAR SYSTEM

For bar constructions



The double retention (internal & external) increases the attachment stability and has superior mechanical properties compared to other attachments. The retention insert maintains constant contact with the abutment, decreasing wear and tear.




COATING

TiN coating for better wear resistance



40° INCLINATION

The Locator Attachment is the lowest vertical attachment. It can counter divergencies up to 40°



MALE PROCESSING PACKAGES



STANDARD
ref. 8519-2

BAR
ref. 8028-2

DIVERGING
ref. 8540-2



Transp. - 2268 grams - ref 8524

Pink - 1361 grams - ref 8527

Blue - 680 grams - ref 8529

STANDARD



Green - 1361 grams - ref 8547

Orange - 450 grams - ref 8915

Red - 230 grams - ref 8548

DIVERGING



Grey - 0 grams - ref 8558

PASSIVE

LOCATOR CORE TOOL: MULTIFUNCTIONAL TOOL



For the full product catalog, contact us via info@proscan.be or via www.proscan.be/webshop